

CAPEX

Carts, Holders, Storage,
Transport & Dispensers



CAPEX

Contents



01

QCars

- 2 Q10 Food Transporter
- 3 Q20 Food Transporter
- 4 Q30 Food Transporter
- 5 Q11 Anaesthetic Cart
- 6 Q11 Procedure Cart
- 7 Q11 Plaster of Paris Cart
- 8 Q11 Emergency Cart
- 9 Q11 Medicine Cart
- 10 Q26 Linen Locker Cart
- 11 Q26 Sterile Supply Cart
- 12 Q11 Nurses'/Doctors' Mobile Work Surface
- 13 Wall-mounted Work Surface
- 14 Acrylic Drawers
- 15 Protective Top & Wire Shelf
- 16 Sub Containers & Divider
- 17 Food & Tea Trays

18

Holders

- 19 Brushes – Wire Mesh Holder
- 20 Gloves – Wire Mesh Holder
- 21 Masks – Wire Mesh Holder
- 22 Hibi Scrub Holder
- 23 D-Germ Bottle Holder
- 24 Caps – Wire Mesh Holder

25

Epoxy-coated Trolleys

- 26 Soiled Linen Theatre Trolley
- 26 Soiled Linen Bag

Epoxy-coated Trolleys Continued

- 27 Make-up/Fixing 2-tier Trolley
- 28 Clean Linen Trolley
- 29 White Cover for Clean Linen Trolley
- 29 General Purpose Heavy Duty Trolley
- 30 General Purpose Trolley
- 31 Anaesthetist Trolley

32

Mobile Drip Stands

- 33 General Mobile Drip Stand
- 34 Theatre Mobile Drip Stand

35

Stainless Steel Trolleys

- 36 Gynae 'Colpo' Trolley
- 37 Small Dressing Trolley
- 38 Medium Dressing Trolley
- 39 Large Dressing Trolley
- 40 Small Surgical Trolley
- 41 Medium Surgical Trolley
- 42 Large Surgical Trolley
- 43 Intubation Trolley
- 44 Admissions Trolley with Holder
- 45 Admissions Trolley
- 46 Single Bowl Stand
- 47 Double Bowl Stand
- 48 Single Patient Wash Trolley
- 49 Double Patient Wash Trolley
- 50 Baby Scale Trolley
- 50 Large Sterile Procedure/Tool Trolley

51

Mayo Table Trolley

- 52 Mayo Table

53

Miscellaneous

- 54 Kick-a-bout
- 54 Stainless Steel Bucket & Lid
- 55 Refrigerator Stand
- 56 Multi-functional LED Timer/Clock

58

Wall-mounted Storage Racks

- 59 Wall-mounted Mop Rack
- 60 Urinal/Bedpan Rack & Drip Tray
- 61 Refuse Bag Holder

62

Mobile Trolley Storage Racks

- 63 Mobile Catheter Storage Rack
- 64 Swab Rack Trolley



Mayo Table

52

65

Urinal Bag Holder Storage

66 Floor-standing Urinal Bag Holder
67 Hanging Urinal Bag Holder

68

Utility Bracket Mountings

69 600 mm Wall-mounted Utility Pillar
69 1200 mm Wall-mounted Utility Pillar

QCARTS

Quality mobile units in a variety
of sizes and functionalities,
complete with acrylic add-ons.

Q10 Food Transporter

CODE: Q10-LOCKER

Q10 transporter/storage locker with one door to a 10-tray compartment, 4 x 125 mm non-marking castors, two swivel castors and two castors with brakes.

ADVANTAGES

- Insulated body manufactured by the thermoforming process in durable ABS polymer with a surface skin of UV & stain-resistant acrylic.
- One door for 10-tray compartment, opening 160° for ease of access – the door is provided with a stainless steel handle.
- Stainless steel push/pull handle located at either end at a convenient height to avoid back strain.
- Moulded wheel base supported by galvanised mild steel frame with an extended edge PVC bumper for protection.
- Includes easy-glide shelves to support the 10 service trays.
- Moulded ABS top has a marine edge to prevent articles sliding off.



FUNCTION

Used for transport of light refreshments from kitchen.



DIMENSIONS

H 1400 mm x W 450 mm x D 600 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Kitchen



Push/pull handle



Easy-glide shelves



Stainless steel handle



Moulded wheel base



Door stopper

Q20 Food Transporter

CODE: Q20-LOCKER

Q20 transporter/storage locker with two doors opening to two 10-tray compartments, 4 x 125 mm non-marking castors, two swivel castors and two castors with brakes.

ADVANTAGES

- Modular, insulated, lightweight but rugged construction.
- Insulated body manufactured by the thermoforming process in durable ABS polymer with a surface skin of UV & stain-resistant acrylic.
- Two doors for two 10-tray compartments, opening 160° for ease of access – each door is provided with a stainless steel handle.
- Stainless steel push/pull handle located at either end at a convenient height to avoid back strain. Moulded wheel base supported by galvanised mild steel frame with an extended edge PVC bumper for protection.
- Includes easy-glide shelves to support the 20 service trays, trays not included.
- Moulded ABS top has a marine edge to prevent articles sliding off.



FUNCTION

The heat retentive material and design make this ideal for the transport of meals/tea to patients or the transport of pharmacy supplies.



DIMENSIONS

H 1400 mm x W 885 mm x D 600 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Kitchen





Q30 Food Transporter

CODE: Q30-LOCKER

Q30 transporter/food locker with three doors, one per 10-tray compartment, 4 x 125 mm non-marking castors, two swivel castors and two castors with brakes. The door of each compartment is fitted with a stainless steel handle to ensure easy opening and closing.

ADVANTAGES

- Modular, insulated, lightweight but rugged construction.
- Insulated body manufactured by the thermoforming process in durable ABS polymer with a surface skin of UV & stain-resistant acrylic.
- Three doors per 10-tray compartment opening 160° for ease of access – each door is provided with a stainless steel handle.
- Stainless steel push/pull handle located at either end at a convenient height to avoid back strain. Moulded wheel base supported by galvanised mild steel frame with an extended edge PVC bumper for protection.
- Includes easy-glide shelves to support the 30 service trays, trays not included.
- Moulded ABS top has a marine edge to prevent articles sliding off.



FUNCTION

The heat retentive material and design make this ideal for the transport of meals/tea to patients or the transport of pharmacy supplies.



DIMENSIONS

H 1400 mm x W 1480 mm x D 600 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Kitchen



Q11 Anaesthetic Cart

CODE: REF.10

Q11 anaesthetic cart with lockable roll front, 4 x 100 mm castors including the following recommended components:

- 1 x protective top
- 1 x AD1 drawer
- 2 x AD2 drawers
- 1 x AD3 drawer
- 24 x sub container dividers
- 12 x long sub containers
- 4 x medium sub containers
- 4 x short sub containers

ADVANTAGES

- Mobile storage (where store rooms are not available).
- Can be pushed wherever it is required.
- Lockable full roller doors.



Lockable

Roll front door

Bumpers

2 castors with brakes



FUNCTION

This trolley is used to store the required medication for anaesthetic in theatre.



DIMENSIONS

H 860 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Theatre Complex



Q11 Procedure Cart

CODE: REF.04

Q11 procedure cart with lockable roll front, 4 x 100 mm castors including the following recommended components:

- 1 x AD1 drawer
- 2 x AD2 drawers
- 1 x AD3 drawer
- 1 x protective top
- 12 x short sub containers

ADVANTAGES

- Mobile storage.
- Can be pushed wherever it is required.
- Lockable full roller doors.



FUNCTION

This trolley is used to store patient medication for issue during ward rounds. Can also be used as a sterile supply cart at ward level. Can also be used in theatre for storing of specialised equipment i.e., ENT, endoscopic equipment, arthroscopy equipment, tracheotomy set.



DIMENSIONS

H 860 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Wards, Out Patient Department



Open cart front with 4 drawers



12 x short sub containers

Lockable
roll front
door

Bucket holder
with stainless
steel bucket



Q11 Plaster of Paris Cart

CODE: REF.11

Q11 plaster of Paris cart with lockable roll front, 4 x 100 mm castors including the following recommended components:

- 1 x AD1 drawer
- 2 x AD2 drawers
- 1 x AD3 drawer
- 1 x protective top
- 12 x medium sub containers
- 24 x sub container dividers
- 1 x bucket holder, stainless steel bucket & lid

ADVANTAGES

- Mobile storage (where store rooms are not available).
- Can be pushed wherever it is required.
- Lockable full roller doors.
- Storage for plaster shear and power tool.



FUNCTION

Used for storage of plaster of Paris and bandages, as well as cast removal equipment. Content can be checked and controlled against the standard stocklist.



DIMENSIONS

H 860 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Theatre Complex



ADVANTAGES
Sufficient space for all requirements. Strong, yet light enough for fast movement.

Q11 Emergency Cart

CODE: REF.03A

Q11 emergency cart with lockable roll front, 4 x 100 mm castors including the following recommended components:

- 1 x AD1 drawer, 2 x AD2 drawers, 1 x AD3 drawer, 1 x protective top
- 12 x large sub containers
- 6 x medium sub containers
- 4 x small sub containers
- 24 x divider sub containers
- 1 x oxygen bottle holder, 1 x CPR board



FUNCTION

Used specifically in emergency situations and is designed in its make up to contain stock, medication and equipment required for resuscitating a patient. Trolley is clearly marked with RED tape and has sufficient space for positioning of the defibrillator.



DIMENSIONS

H 860 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

ER, Trauma Unit



Protective top



Detachable CPR board

RECOMMENDED STOCKLIST

TOP SHELVE:

- Cold Light Source
- Fibre Optic Bronchoscope
- Channel Cleaning Brush

DRAW NO. 1:

- Laryngoscope Handle
- Laryngoscope Blades No. 1
- Laryngoscope Blades No. 2

- Laryngoscope Blades No. 3
- Laryngoscope Blades No. 4
- Laryngoscope Blades Miller No. 2
- Xylocaine Spray with Nozzle
- Forceps Artery Crile Straight 15 cm
- Forceps Introducing McGills
- Adult 20 cm
- Primapore Small
- Silk Tape Small



12 x short sub containers



Lockable roll front door

Q11 Medicine Cart on Height Extension

CODE: REF.15

Q11 medicine cart with lockable roll front, 4 x 100 mm castors, 120 mm epoxy-coated extended base including the following recommended components:

- 1 x AD1 drawer
- 2 x AD2 drawers
- 1 x AD3 drawer
- 1 x protective top
- 12 x short sub containers

ADVANTAGES

- Mobile storage (where store rooms are not available).
- Can be pushed wherever it is required.
- Lockable full roller doors.



FUNCTION

Used to store patient medication for issue during ward rounds.



DIMENSIONS

H 980 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Wards

- Transpore Small
- Syringe 20 cc
- Syringe 5 cc
- Syringe 2 cc
- Needles 23 g Pink x 4
- Remicaine Jelly Tube x 1
- KY Jelly Tube x 1
- Bite Block Blachey x 1
- Bite Block Santek x 5

- Bite Block Rubber x 1
- Tube Aspiration Luki
- Crico-Cricothyroidotomy Kit
- Laryngotracheal Anaesthesia Kit

DRAW NO. 2:

- Introducer Stylet
- Introducer Stylet Extra Long

- Introducer Stylet Paediatric
- Introducer Stylet Rüsh
- Introducer
- Stylet Satin Slip x 2
- Airway Nasopharyngeal Size 6.0
- Airway Nasopharyngeal Size 7.0
- Airway LMA Fasträch Size 3

- Airway LMA Fasträch Size 4
- Airway LMA Fasträch Size 5
- Stabiliser Fasträch x 1
- Laryngeal Mask Size 2
- Laryngeal Mask Size 2½
- Laryngeal Mask Size 3
- Laryngeal Mask Size 4
- Laryngeal Mask Size 5

DRAW NO. 3:

- Laryngoscope McCoy
- Blade No. 1
- Blade No. 2
- Blade No. 3
- Blade No. 4



Q26 Linen Locker Cart

CODE: REF.02A

Q26 linen locker with lockable roll fronts,
4 x 125 mm castors including:

- 5 x shelves of stainless steel wire.

ADVANTAGES

- Mobile storage (where store rooms are not available).
- Can be pushed wherever it is required.
- Control of linen stock.
- Lockable full roller doors.



FUNCTION

This trolley is used to transport both ward and theatre linen from the laundry to the point of requirement. Shelves are stainless steel and allow air circulation for linen to remain fresh.



DIMENSIONS

H 1750 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Wards, Out Patient Department



Lockable roll front doors



Push/pull handles



Q26 Sterile Supply Cart

CODE: REF.06

Q26 sterile supply with lockable roll fronts, 4 x 125 mm castors including the following recommended components:

- 1 x AD1 drawer
- 5 x AD2 drawers
- 2 x AD3 drawers
- 10 x long sub containers
- 6 x medium sub containers
- 10 x short sub containers
- 20 x sub container dividers

ADVANTAGES

- Mobile storage (where store rooms are not available).
- Can be pushed wherever it is required.
- Cost control and cost recovery is easy as a content list can be placed inside the unit.



FUNCTION

Used at ward and theatre level, specifically for the storage of sterile supplies or other ward stock.



DIMENSIONS

H 1750 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Theatre Complex, Out Patient Department



Q11 Nurses'/Doctors' Mobile Work Surface

CODE: REF.18

Q11 cart without lockable roll front, on 4 x 100 mm castors, 120 mm extended base, including the following recommended components:

- 1 x AD1 drawer
- 1 x AD2 drawer
- 1 x AD3 drawer
- 12 x short sub containers

ADVANTAGES

- Mobile writing surface for circulating staff or doctors.
- Storage drawers for pathology forms, general needs during operating lists, e.g., linen savers, operating table covers.
- Easy to handle.



FUNCTION

Used as writing surface in operating room. Also, where dry walls are not suitable for wall-mounted nurse'/doctors' writing surfaces.



DIMENSIONS

H 980 mm x W 610 mm x D 610 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Drawers, sub containers, top and bottom roller doors, protective tops, castors and bumpers are available on request.



DEPARTMENT

Theatre Complex



R-2 Wall-mounted Nurses'/Doctors' Work Surface

CODE: 2009-27

Moulded C-frame including:

- 1 x AD1 drawer
- 1 x AD2 drawer
- 1 x AD3 drawer
- 12 x short sub containers

ADVANTAGES

- Storage drawers for pathology forms, general needs during operating lists, e.g. linen savers, operating table covers.



FUNCTION

Used as writing surface in operating room.



DIMENSIONS

W 610 mm x D 580 mm HEIGHT FROM FLOOR LEVEL: 930 mm




CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



#1

Acrylic Drawer #1 Small

CODE: 260-200-AD1

Transparent, vacuum formed drawer manufactured from PETG material (polyethylene copolymer). Sterilisable and recyclable.



FUNCTION

Used as storage compartments in combination with the QCart storage and transport system.



DIMENSIONS

L 540 mm x W 510 mm x H 40 mm



CLEANING

To be cleaned with a solution of warm water with neutral soap, rinsed with water using a soft sponge. To be preserved under normal operating conditions (under 300°C).



#2

Acrylic Drawer #2 Medium

CODE: 260-200-AD2


DIMENSIONS

L 540 mm x W 510 mm x H 100 mm




#3

Acrylic Drawers #3 Large

CODE: 260-200-AD3


DIMENSIONS

L 540 mm x W 510 mm x H 150 mm



Protective Top

CODE: 260-100-PTP

Transparent, vacuum formed protective top, with raised-edge design, manufactured from PETG material (polyethylene copolymer). Sterilisable and recyclable.



FUNCTION

Used as cover to protect QCart top from damage. Prevents stationery placed on QCart top from rolling off.



DIMENSIONS

L 615 mm x W 620 mm



CLEANING

To be cleaned with a solution of warm water with neutral soap, rinsed with water using a soft sponge. To be preserved under normal operating conditions (under 300°C).

Wire Shelf

CODE: 260-16-WSH_A

Shelf wire manufactured from stainless steel (grade 304) frame with electropolish finish. Allow air circulation for linen to remain fresh.



FUNCTION

Used in combination with Q26 cart (Ref.02A) for storage of linen, transported from laundry to the point of requirement.



DIMENSIONS

L 538 mm x W 506 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Kitchen/Cafeteria

Sub Container Short

CODE: 260-300-SCS




DIMENSIONS

L 225 mm x W 75 mm x H 75mm



Sub Container Medium

CODE: 260-300-SCM



DIMENSIONS

L 225 mm x W 150 mm x H 75 mm



Sub Container Long

CODE: 260-300-SCL



DIMENSIONS

L 450mm x W 75 mm x H 75 mm



Sub Container Divider

CODE: 260-300-075



DIMENSIONS

L 66 mm x W 45 mm



Sub Containers

Transparent containers manufactured from acrylic material.



FUNCTION

Used to perfectly store smaller objects/sterile supplies.



CLEANING

To be cleaned with a solution of warm water with neutral soap, rinsed with water using a soft sponge. To be preserved under normal operating conditions (under 300°C).

Sub Container Divider

Transparent dividers manufactured from acrylic material.



FUNCTION

Used to divide sub containers to separate required medication/sterile supplies.



CLEANING

To be cleaned with a solution of warm water with neutral soap, rinsed with water using a soft sponge. To be preserved under normal operating conditions (under 300°C).

Food Tray

2012-260

Strong
and easy to
clean



Crockery & cutlery not included

Tea Tray

2012-262



Crockery & cutlery not included

Food & Tea Trays

ABS with a slip-resistant finish on working face and four dimples on reverse side of tray. Allows a lateral tilt of $\pm 15^\circ$ with full cup and saucer before slipping.

ADVANTAGES

- Dimples on reverse side of tray give stability and allows airflow when stacked.
- Takes up minimal space when not in use.



FUNCTION

Used as service tray for food, tea/coffee.



DIMENSIONS

Food tray: L 540 mm x W 386 mm
Tea tray: L 296 mm x W 215 mm



CLEANING

Use gentle detergents for cleaning. Do not use steel wool/scouring pads for cleaning.



DEPARTMENT

Kitchen/Cafeteria

18

CAPEX

HOLDERS

A wide variety of stainless steel holders for various uses.

Innovating
since 1988.



Brushes – Stainless Steel Wire Mesh Holder

CODE: 29-00

Brush holder constructed with 5 mm stainless steel (grade 304) round bar, with electropolish finish.

ADVANTAGES

- Stainless steel, therefore durable.
- Removable.
- Easy to clean.
- Visual control over stock.
- Can be washed in instrument washer.
- Autoclavable.



FUNCTION

Used as storage for sterile brushes at scrub up area.



DIMENSIONS

H 615 mm x W 150 mm x D 145 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theater Complex



Scrub
brushes not
included



Gloves – Wire Mesh Stainless Steel Holder

CODE: 32-00

Wire mesh glove holder consisting of 5 mm stainless steel (grade 304) round bar with electropolish finish.

ADVANTAGES

- Stainless steel.
- Removable.
- Easy to clean.
- Visual control over stock.
- Can be washed in instrument washer.
- Autoclavable.



FUNCTION

Supply for daily use of surgical gloves (different sizes).



DIMENSIONS

H 676 mm x W 160 mm x D 120 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Masks – Wire Mesh Stainless Steel Holder

CODE: 30-00

Wire mesh mask holder consisting of 5 mm stainless steel (grade 304) round bar with electropolish finish.

ADVANTAGES

- Stainless steel.
- Removable.
- Easy to clean.
- Visual control over stock.
- Can be washed in instrument washer.
- Autoclavable.



FUNCTION

Supply for daily use of surgical face masks.



DIMENSIONS

H 208 mm x W 124 mm x D 132 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Stainless Steel Hibi Scrub MK2 Holder

CODE: 341-00

Stainless steel (grade 304) holder with electropolished finish. Designed to accommodate a round 500 ml bottle.

ADVANTAGES

- Stainless steel, therefore durable.
- Removable.
- Easy to clean.
- Can be washed in instrument washer.
- Autoclavable.



FUNCTION

Used at surgical scrub basins as hands-free soap dispenser.



DIMENSIONS

H 245.5 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex, Gastro
Enterology, Out Patient Department



Stainless Steel D-Germ Bottle Holder

CODE: 343-00

Holder for alcohol/surface cleaning bottle, manufactured from stainless steel (grade 304) with electropolish finish.

ADVANTAGES

- Stainless steel, therefore durable.
- Removable.
- Easy to clean.
- Can be washed in instrument washer.
- Autoclavable.



FUNCTION

Used in various clinical areas as bottle holder e.g. patient rooms.



DIMENSIONS

H 271 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards, Critical Care, High Care, Neonatal, ER, Theatre Complex, Laundry, Pharmacy

Bottle not
included





Caps – Wire Mesh Stainless Steel Holder

CODE: 31-00

Wire mesh cap holder consisting of 5 mm stainless steel (grade 304) round bar with electropolish finish.

ADVANTAGES

- Stainless steel.
- Removable.
- Easy to clean.
- Visual control over stock.
- Can be washed in instrument washer.
- Autoclavable.



FUNCTION

Supply for daily use of surgical caps.



DIMENSIONS

H 275 mm x W 205 mm x D 205 mm



CLEANING

To be cleaned with a disinfectant agent,
NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



CAPEX

25

EPOXY- COATED TROLLEYS

A wide variety of sizes and styles for
an extensive range of functions.

Always
client-
centred.



Soiled Linen Theatre Trolley Excluding Bag

CODE: 25-00

Soiled linen trolley consisting of epoxy-coated u-shaped steel frame (colour: motor craft grey) with four hooks (for bag) and 4 x 100 mm swivel castors.



FUNCTION

Receiver trolley for collection of soiled linen.



DIMENSIONS

H 1050 mm x W 700 mm x D 370 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Bag attachment hooks



Brass eyelets for hooks

Soiled Linen Bag

CODE: 25-10



White 100% polyester bag including 4 x brass eyelets, to be hooked onto soiled linen trolley (Ref.25-00).



FUNCTION

Bag available with soiled linen receiver trolley.



DIMENSIONS

H 820 mm x W 390 mm x L 610 mm



CLEANING

Machine washed at 60°C. Tumble dry on low heat.



Make-up/Fixing 2-tier Trolley

CODE: 27-01

Trolley consisting of mild steel frame, white epoxy coated, including 2-tier plastic contiki-slatted shelves, including 4 x 125 mm swivel castors, 4 x 95 mm grey bumpers.



FUNCTION

Transportation of clean linen to the patient's bed during linen changes.



DIMENSIONS

H 990 mm x W 915 mm x L 438 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



Clean Linen Trolley

CODE: 26-10

Trolley for clean linen – 5-tier with epoxy-coated steel frame (colour: motor craft grey), 4 x 125 mm castor wheels and 10 x plastic contiki-slatted shelves. Safety – maximum load-bearing ± 15 kg per shelf.

ADVANTAGES

- Very mobile and easy to clean.
- Accompanied by zipped canvas cover, labelled appropriately.
- 10-removable slats for easy cleaning.



FUNCTION

Extremely versatile, ideally suited to transport clean linen between laundry and wards.



DIMENSIONS

H 1851 mm x W 915 mm x D 645 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



Wheel size can be adapted to accommodate uneven terrain (Ref. 26-12)



White Cover for Clean Linen Trolley

CODE: 26-11

White 100% polyester cover, 1 x blue zip each side, including inscription 'Clean Linen'.



FUNCTION

Use to cover clean linen trolley during the transportation from laundry to nursing units.



DIMENSIONS

H 1615 mm x W 910 mm x D 450 mm



CLEANING

Machine washed at 60°C. Tumble dry on low heat.



DEPARTMENT

Wards

General Purpose Heavy Duty Trolley

CODE: 26-12

Trolley for clean linen – 5-tier with epoxy-coated steel frame (colour: motor craft grey), 4 x 200 mm wheelchair castors and 10 x plastic contiki-slatted shelves.



FUNCTION

Extremely versatile, ideally suited to transport clean linen between laundry and wards over uneven terrain.



DIMENSIONS

H 1566 mm x W 985 mm x D 645 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



Extremely
versatile

General Purpose Trolley

CODE: 26-00

General purpose trolley – 5-tier with epoxy-coated steel frame (colour: motor craft grey), 4 x 125 mm swivel castors, 4 x plastic bumper wheels and 10 x plastic contiki-slatted shelves.



FUNCTION

Use as storage space. To transport sets from sterile store to sterile set up.



DIMENSIONS

H 1566 mm x W 985 mm x D 645 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex

10-removable slats for easy cleaning



Anaesthetist Trolley

CODE: 02-01

Epoxy-coated steel frame with following components: top, drawer, 4 x 75 mm castors with circlip stems (two with brakes and two without brakes).

ADVANTAGES

- Mobile with two castors with brakes.
- Drawer for personal belongings.
- Shelf for briefcase.
- Right or left handed use.



FUNCTION

Spacious writing surface for anaesthetist in operating room.



DIMENSIONS

H 600mm x W 950mm x D 450 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



32

CAPEX

MOBILE DRIP STANDS

Safe, durable and versatile.

Always
prioritising
patient
safety.





Adjustable height



4-hook bag holder



Counter weight for stability



Durable 75 mm twin wheel castor

**MAXIMUM
LOAD CAPACITY:
10 KG**

General Mobile Drip Stand

CODE: 14-135

Mobile drip stand consisting of a 5-star counter-weighted glass-filled nylon base with 5 x 75 mm castors with a 25 mm diameter pillar and a stainless steel (grade 304) top pole 19 mm in diameter to which 4 x stainless steel (grade 304) drip hooks are attached.

ADVANTAGES

- Stainless steel (grade 304) tubing is used instead of chrome-plated tubing.
- The 4-hook bag holder can be adjusted lower and/or higher (longer regulating range).
- Whilst releasing the knob for lowering the bag holder, especially when loaded, the bag holder cannot accidentally drop on the operator's hand as it will stop 100 mm from the release knob.
- For increasing the stability, especially when pumps are attached to the stand, the weight of the base is increased with the addition of the counter weight (total weight of 8 kg).
- An aesthetically pleasing, easy to clean, glass-filled nylon 5-star base is used.
- A durable 75 mm twin wheel castor is being used to improve the mobility of the stand to prevent a topple over when it reaches obstacles. i.e. carpet strips.



FUNCTION

Support for intravenous vaculitres.



DIMENSIONS

MAXIMUM HEIGHT: 2170 mm
MINIMUM HEIGHT: 1250 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



OPTIONAL COMPONENTS

Basket on pillar (code: IHB-04-05A)



DEPARTMENT

Wards, ER

**MAXIMUM
LOAD CAPACITY:
10 KG**

Theatre Mobile Drip Stand

CODE: 2013-14-135

Mobile drip stand consisting of a 5-star counter-weighted glass-filled nylon base with 5 x 75 mm castors with a 25 mm diameter pillar and a stainless steel (grade 304) top pole 19 mm. 500 mm longer than the drip stand used in wards, to which 4 x stainless steel (grade 304) infusion drip hooks are attached.

ADVANTAGES

- The drip stand top pole can be extended 500 mm higher when a patient is resuscitated.
- In case of an emergency, the hooks are more accessible and user friendly when IV bags are attached.
- Stainless steel (grade 304) tubing is used instead of chrome-plated tubing.
- Easy to clean stainless steel 4-hook bottle holder instead of plastic.
- The 4-hook bag holder can be adjusted lower and/or higher (longer regulating range).
- Whilst releasing the knob for lowering the bag holder, especially when loaded, the bag holder cannot accidentally drop on the operator's hand as it will stop 100 mm from the release knob.
- For increasing the stability, especially when pumps are attached to the stand, the weight of the base is increased with the addition of the counter weight (total weight of 8 kg)
- Aesthetically pleasing, easy to clean, glass-filled nylon 5-star base is used.
- A durable 75 mm twin wheel castor is being used to improve the mobility of the stand to prevent a topple over when it reaches obstacles i.e. carpet strips.



FUNCTION

Support for intravenous vaculitres.



DIMENSIONS

MAXIMUM HEIGHT: 2670 mm

MINIMUM HEIGHT: 1250 mm



CLEANING

To be cleaned with a disinfectant agent,
NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex

CAPEX



35

STAINLESS STEEL TROLLEYS

Choose from a variety of sizes to cater
to a wide range of functions.

Premium
quality
medical
solutions.

Gynae 'Colpo' Trolley

CODE: 155-00

Stainless steel (grade 304) frame and top,
4 x 75 mm twin wheel swivel castors.



FUNCTION

Used during surgical
gynaecological lithotomy and
urology procedures as sterile
work surface to surgeon.



DIMENSIONS

H 738 mm x W 827 mm x D 430 mm



CLEANING

To be cleaned with a disinfectant agent,
NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Swivel castors with brakes



Pads to avoid vibration noise

Small Stainless Steel Dressing Trolley

CODE: 15-DS

Dressing trolley consisting of stainless steel (grade 304) frame including a stainless steel rim on top shelf on three sides, with two removable stainless steel (grade 304) shelves, and 4 x 75 mm twin wheel castors with circlip stems.



FUNCTION

Used for patient bedside dressing or treatment.



DIMENSIONS

H 895 mm x W 435 mm x D 415 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Ideally for doctors' rooms, small theatres and clinics.



Medium Stainless Steel Dressing Trolley

CODE: 16-DM

Dressing trolley consisting of stainless steel (grade 304) frame including a stainless steel rim on top shelf on three sides, with two removable stainless steel (grade 304) shelves, and 4 x 75 mm twin wheel castors with circlip stems.



FUNCTION

Used for patient bedside dressing or treatment.



DIMENSIONS

H 895 mm x W 565 mm x D 415 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Ideally for doctors' rooms, small theatres and clinics.



Stainless steel (grade 304)



Swivel castors with bumpers



Pads to avoid vibration noise.



Large Stainless Steel Dressing Trolley

CODE: 17-DL

Dressing trolley consisting of stainless steel frame (grade 304) including a stainless steel rim on top shelf on three sides, with two removable stainless steel (grade 304) shelves, and 4 x 75 mm twin wheel castors with circlip stems.



FUNCTION

Used for patient bedside dressing or treatment.



DIMENSIONS

H 895 mm x W 900 mm x D 435 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Ideally for doctors' rooms, small theatres and clinics.



Small Stainless Steel Surgical Trolley

CODE: 18-IS

Surgical trolley consisting of stainless steel (grade 304) frame, with two removable stainless steel (grade 304) shelves, and 4 x 75 mm twin wheel castors with circlip stems.



FUNCTION

Serves as work surface during sterile surgical procedures in operating rooms.



DIMENSIONS

H 862 mm x W 450 mm x D 450 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Pads to avoid vibration noise.



Swivel castors with bumpers.





MAXIMUM
WEIGHT
CAPACITY:
80 KG

Medium Stainless Steel Surgical Trolley

CODE: 19-SM

Surgical trolley consisting of stainless steel (grade 304) frame with two stainless steel (grade 304) shelves and 4 x 75 mm twin wheel castors with circlip stems.

ADVANTAGES

- Improved double wheel castors.
- Improved stainless steel polish finish throughout.
- Closer fabrication tolerances on sheet metal & framework.
- Bump-on pads to avoid vibration noise.



FUNCTION

Only used during sterile surgical procedures in operating room.



DIMENSIONS

H 860 mm x W 600 mm x D 450 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex/Gastro Enterology

MAXIMUM
WEIGHT
CAPACITY:
100 KG



Large Stainless Steel Surgical Trolley

CODE: 20-SL

Surgical trolley consisting of stainless steel (grade 304) frame with two stainless steel (grade 304) shelves and 4 x 75 mm twin wheel castors with circlip stems.

ADVANTAGES

- Improved double wheel castors.
- Improved stainless steel polish finish throughout.
- Closer fabrication tolerances on sheet metal & framework.
- Bump-on pads to avoid vibration noise.



FUNCTION

Only used during sterile surgical procedures in operating room.



DIMENSIONS

H 860 mm x W 915 mm x D 450 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex

**MAXIMUM
WEIGHT
CAPACITY:
80 KG**

Intubation Trolley

CODE: 21-IS

Trolley consisting of stainless steel (grade 304) frame with two stainless steel (grade 304) shelves, two medium acrylic drawers, 4 x 75 mm twin wheel castors with circlip stems, two rails fitted to both shelves and two hooks for anaesthesia screen and a resuscitation bag.

ADVANTAGES

- Improved twin wheel castors.
- Improved stainless steel polish finish throughout.
- Closer fabrication tolerances on sheet metal & framework.
- Bump-on pads to avoid vibration noise.



FUNCTION

Only used for storage of intubation supplies.



DIMENSIONS

H 903 mm x W 600mm x D 510 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Slide out work surface



Rail for anaesthesia screen



Hook for resuscitation bag



MAXIMUM
WEIGHT
CAPACITY:
80 KG

Admissions Trolley with D-Germ Bottle Holder & Sharps Container Holder

CODE: 24-AD

Trolley consisting of stainless steel rim on top and bottom shelves on three sides, with two removable stainless steel (grade 304) shelves, two medium acrylic drawers, pull out shelf and 4 x 75 mm twin wheel castors with circlip stems.



FUNCTION

Used by nursing staff during the admission of patients to the ward.



DIMENSIONS

H 231 mm x W 231 mm x D 161 mm
D-Germ Bottle Holder: H 143 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



D-Germ Bottle Holder &
Sharps Container Holder



Admissions Trolley

CODE: 22-AD

Trolley consisting of stainless steel (grade 304) frame, with two removable stainless steel (grade 304) shelves, two medium acrylic drawers, pull out shelf and 4 x 75 mm twin wheel castors with circlip stems.



FUNCTION

Used by nursing staff during the admission of patients to the ward.



DIMENSIONS

H 860 mm x W 450 mm x D 600 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



Single Stainless Steel Bowl Stand without Bowl

CODE: 234-00

Stainless steel (grade 304) frame on 4 x 75 mm twin wheel castors with circlip stems, including four bumper wheels. Bowl not included.



FUNCTION

Use during surgical procedures for hand rinsing.



DIMENSIONS

H 923 mm x W440 mm x D 460 mm
Diameter of ring: 300 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Bowl stand with bowl



Excludes
stainless steel
bowls

Double Stainless Steel Bowl Stand without Bowl

CODE: 234-01

Stainless steel (grade 304) frame on 4 x 75 mm twin wheel castors with circlip stems, including four bumper wheels.



FUNCTION

Use during surgical procedures for hand rinsing.



DIMENSIONS

H 923 mm x W 625 mm x L 460 mm
Diameter of ring: 300 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex





Single Patient Wash Trolley without Bowl

CODE: 238-01

Trolley consisting of stainless steel (grade 304) frame, with single removable contiki-slatted shelf, colour white, including 4 x 75 mm twin wheel castors with circlip stems, 4 x 95 mm grey bumpers. Excluding wash bowl.



FUNCTION

Used for washing/bathing of patients.



DIMENSIONS

H 931 mm x W 509 mm x D 460 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



Bowl stand with bowl



Double Patient Wash Trolley without Bowls

CODE: 239-00

Trolley consisting of stainless steel (grade 304) frame, with two removable contiki-slatted shelves, colour white, including 4 x 75 mm twin wheel castors with circlip stems, 4 x 95 mm grey bumpers. Excluding wash bowls.



FUNCTION

Used for washing/bathing of patients.



DIMENSIONS

H 928 mm x W 922 mm x D 460 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards

Baby Scale Trolley

CODE: 2012-14

Trolley consisting of stainless steel frame including one removable stainless steel shelf, including 4 x 75 mm twin wheel castors with circlip stems, 4 x 95 mm grey bumpers.



FUNCTION

Storage space for baby scale (weighing of babies at a practical height).



DIMENSIONS

H 870 mm x W 662 mm x D 400 mm




CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Paediatric Unit/ER



Large Sterile Procedure/Tool Trolley

CODE: 23-SP

Stainless steel (grade 304) frame on 4 x 75 mm twin wheel castors with circlip stems, including four bumper wheels.



FUNCTION

Use during surgical procedures for hand rinsing. Bowl not included.



DIMENSIONS

H 860 mm x W 915 mm x D 450 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



CAPEX

51

MAYO TABLE TROLLEYS

Sterile, durable and suitable
for any operating room.

Driven by
performance.



Mayo Table

CODE: 325-00

Stainless steel (grade 304) trolley consisting of the following components: stainless steel (grade 304) frame that fits over base of operating table, 4 x 95 mm bumper wheels and 4 x 75 mm twin wheel castors with circlip stems (two with brakes and two without brakes).



FUNCTION

Used during surgical procedures as sterile working surface. Height adjustable.



DIMENSIONS

W 527 mm x D 372 mm BASE: W 700 mm x D 406 mm
MINIMUM HEIGHT: 1030 mm (floor to tray)
MAXIMUM HEIGHT: 1300 mm (floor to tray)



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



MISCELLANEOUS

From stainless steel buckets
to refrigerator stands.

Always
team-
orientated.

Kick-a-Bout

CODE: 2013-228

The kick-a-bout is easily manoeuvrable and functional, consisting of a rigid and durable white urethane ring, fitted on 4 x 50 mm twin wheel castors with circlip stems. Can be easily pushed from one place to another with just a 'kick'.



FUNCTION

Secure holder for bucket for the receiving of medical waste.



DIMENSIONS

OUTER DIAMETER: 420 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Stainless Steel Bucket & Lid

CODE: SSBL-01

Stainless steel 12-litre bucket with lid.
Accessory to 2013-228.



DIMENSIONS

DIAMETER: 305 mm x H 280 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex





Refrigerator Stand

CODE: 47-00

Stand consisting of mild steel frame, white epoxy coated, including four levelling feet to lend a steady position to the refrigerator.



FUNCTION

Stand for medicine refrigerator.



DIMENSIONS

H 630mm x W 474 mm x D 434 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



Multi-functional LED Timer/Clock

CODE: CTS6D4N-BR-00

Elapsed timer solution includes a LED digital clock, used in combination with a start/stop control panel, to start or stop the counting function; beeper, and hour and minute reset buttons. During emergencies, hospital personnel may activate the timer to start the count with the push of a button on the start/stop control panel.

ADVANTAGES

- Count up (elapsed time)/count down (time remaining)
- Standard clock display
- Start/stop button starts or stops timing session
- 6-digit 100 mm display format
- LED display colour: red
- DC 12V/2A power supply
- Constructed from a matt aluminium alloy case and matt black lens
- Infrared remote control (excluding 2 x AAA batteries)



FUNCTION

Used to allow nurses, doctors and emergency staff to perform life-saving treatments during emergency situations, in a timely manner.



DIMENSIONS

H 160 mm x L 640 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Perfect solution for areas of a healthcare facility such as:

- Operating Rooms, Emergency Rooms, Post-operative Area, Cardiac Areas, Intensive Care Units

Can also be used for time-sensitive events such as:

- Educational Testing, Seminars, Manufacturing Processes, Lab Testing, Recreational Activities





58

CAPEX

WALL- MOUNTED STORAGE RACKS

For all your storage needs.

Proud to be
a level 4 BEE
contributor.





Stainless Steel Wall-mounted Mop Rack

CODE: 168-00

Consisting of stainless steel (grade 304) tubing with 5 x 40 mm stainless steel (grade 304) 6 mm round bar hooks, with stainless steel polish finish. Wall mounted via stainless steel plate.



FUNCTION

Used as a mop rack in sorting/sluice in CSSD area.



DIMENSIONS

L 300 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex, Gastro Enterology, Out Patient Department

Off the floor
storage



Stainless Steel Urinal/Bedpan Rack Including Drip Tray

CODE: 65-00

Consisting of stainless steel (grade 304) bedpan and urinal rack including two removable stainless steel (grade 304) drip trays.



FUNCTION

Storage for six bedpans and seven urinals.



DIMENSIONS

H 818 mm x W 820 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex, Wards, ER, Trauma Unit, Critical Care, Gastro Enterology, Out Patient Department



Refuse Bag Holder

CODE: 275-00

Refuse bag holder consisting of epoxy-coated steel frame with support to keep refuse plastic bag firmly in place.



FUNCTION

To be used at washbasins for general waste e.g. used paper towels. Off the floor for cleaning. Easy and accessible for use and changing of bag.



DIMENSIONS

W 435 mm x D 247 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex, Gastro Enterology, Out Patient Department

62

CAPEX

MOBILE TROLLEY STORAGE RACKS

Quality mobile storage for
your theatre complex.

Set apart by
our innovative
solutions.





Mobile Catheter Storage Rack

CODE: 35-01

Mobile stainless steel (grade 304) frame, containing 16 catheter hooks, 4 x 95 mm bumper wheels and 4 x 75 mm twin wheel castors with circlip stems.

ADVANTAGES

- Quick access and retrieval.
- Visual inventory.
- Easy identification.



FUNCTION

Providing storage solutions for catheters.



DIMENSIONS

H 1200 mm x W 1525 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex



Swab Rack Trolley

CODE: 59-04

Trolley consisting of stainless steel (grade 304) frame including removable drip tray, including 4 x 75 mm twin wheel castors with circlip stems, 4 x 95 mm grey bumpers. Excluding white board.



FUNCTION

Management of swab count in operating theatre.



DIMENSIONS

H 964 mm x W 867 mm X D 477 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.




DEPARTMENT

Theatre Complex

URINAL BAG HOLDER STORAGE

For the comfort and discretion
of your patients.

Commitment
to continuous
improvement.



Floor Standing Urinal Bag Holder

CODE: 179-00

Consisting of stainless steel (grade 304), 6 mm round bar with electropolish finish.



FUNCTION

Used as standing holder for urinal bag.



DIMENSIONS

H 510 mm x W 194 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex, Wards





Hanging Urinal Bag Holder

CODE: 179-01

Consisting of stainless steel (grade 304) bracket, 6 mm round bar with electropolish finish.



FUNCTION

Used as hanging holder for urinal bag.



DIMENSIONS

H 106 mm x W 151 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Theatre Complex, Wards

68

CAPEX

UTILITY BRACKET MOUNTINGS

Conveniently available in two lengths.

A legacy of
innovation.



600 mm Wall-mounted Utility Pillar

CODE: : IHB-19

A wall-mounted utility pillar manufactured from stainless steel (grade 304).



FUNCTION

Mounted in sterile supply store for the storage of infusion pumps.



DIMENSIONS

L 600 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



1200 mm Wall-mounted Utility Pillar

CODE: IHB-23

A wall-mounted utility pillar manufactured from stainless steel (grade 304).



FUNCTION

Mounted in sterile supply store for the storage of infusion pumps.



DIMENSIONS

L 1200 mm




CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards



T: +27 (0)21 851 8484

E: sales@medicalinnovations.co.za

W: medicalinnovations.co.za

7 Cyclonite St, The Interchange
Somerset West, 7130
Western Cape, South Africa

The information and dates supplied in this catalogue, are based on the information believed to be reliable, however are not binding and inaccuracies may occur. Medical Innovations, reserves the right to change specification and design at any time without notice.

© 2017 by South Africa. All rights reserved under the Copyright Act of 1987 with the amendments. Reproduction of any of the content, full or partial, is allowed with prior written approval of South Africa. No reasonable requests will be refused provided the condition as laid out are met.