

MEDICAL DEVICES & CONSUMABLES

Masks, Medical Devices, Gloves,
Protective Clothing and Cleaning & Disinfecting



MEDICAL DEVICES
& CONSUMABLES

Contents



1 Masks

- 1 KN95 Protective Respirators FFP2
- 1 Surgical Face Mask 3-ply with Ear Loops

2 Medical Devices

- 2 ThermoFlash Infrared Handheld Thermometer
- 3 Surgical Clipper Rechargeable

4 Gloves

- 4 MICARE Nitrile Examination Gloves
- 7 Nitrile Examination Gloves

8 Protective Clothing

- 8 Protective Coverall
- 9 Isolation Gown

10 Linen Savers

- 10 Linen Savers

12 Cleaning & Disinfecting

- 12 Surgical Scrub Brush with Nail Cleaner Sterile

13 Dispenser

- 13 Manual Alcohol/Soap Dispenser including Drip Tray & Pump
- 14 Manual Alcohol/Soap Dispenser Only
- 15 Manual Alcohol/Soap Dispenser with Refill Containers
- 16 Alcohol/soap Touchless Dispenser
- 17 Liquid Hand Sanitiser 800 ml Sachet
- 17 Dispenser Pump
- 17 Dispenser Refill Container
- 17 Dispenser Drip Tray

KN95 Protective Respirators FFP2

CODE: MASK-DFMS-KN95

Multiple protection, more than 95% filter to block harmful substances, truly clean breathing. For total peace of mind and safe breathing.

ADVANTAGES

- Anti-bacterial
- Anti-dust
- Anti-fume
- Anti-exhaust
- 3D shape design
- Quadruple filtering
- Multi-layered protection
- High quality
- Super soft non-woven fabric
- Ultrasonic spot-welding technology



QUANTITY
20pcs/pack



Surgical Face Mask 3-ply with Ear Loops

CODE: MASK-MED0005

Used by healthcare professionals to protect themselves and others from the spread of airborne infectious diseases and bodily fluids. Masks are generally designed to completely cover both the mouth and nose.

ADVANTAGES

Masks prevent large droplets of bodily fluids that may contain viruses or germs from escaping from the nose and mouth such as from someone sneezing and coughing.



QUANTITY
50pcs/pack





ThermoFlash Infrared Handheld Thermometer

CODE: MI-THER-LX26E

ThermoFlash LX-26 Evolution is a non-contact thermometer. It is equipped with advanced MicroSecondFlash™ technology, a faster auto-calibration system that reduces error margins. Direct keys without setting functions facilitate use. Fast, it measures the body temperature at a distance of 5 cm. Very hygienic and non-invasive, it detects the heat emitted by the temporal artery (blood flow). By combining technical expertise, precision and simplicity, the LX-26 Evolution is ideal for the whole family's temperature taking.

ADVANTAGES

- Contactless thermometer
- Measure forehead temperature
- Temperature taking in all circumstances, even in low-light environment (back lighted LCD screen)
- Sound alarm will alert you in case of excessive temperature



FUNCTION

Simply direct the ThermoFlash Infrared Thermometer towards the patient's forehead and read the temperature displayed on the LCD screen. As there is no contact, it is perfectly hygienic.



DIMENSIONS

H 175 mm x W 110 mm x D 48 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Wards

Surgical Clipper Rechargeable

CODE: SC001

The world's first surgical clipper with a twistable shaving head.

ADVANTAGES

For Patients

- Minimal risk of micro-injuries, infection, skin irritations
- Low noise level
- Quick and simple procedure, particularly for children
- Eliminates painful cuts, particularly in the genital area
- Better adhesion of wound dressing materials

For Clinicians

- Efficient skin preparation for surgery: time saving, quick and easy
- Safe removal of body hair without cuts and skin lesions
- A clean and clear incision site, facilitating adhesion of dressings to the body
- Significant reduction in infection rates

BENEFITS:

- Efficient body hair removal in both wet and dry surfaces
- The right blade for the right type of clipping: a universal blade and a neuro blade
- High efficiency and thorough shaving
- A high battery performance and short charging time

FUNCTION

Twistable head: safe operation of the clippers in hard to reach areas of the body. The right blade for the right type of clipping:

- Universal blade and a neuro blade
- Non-slip grip - no risk of injury from slipping
- Easy to use: insert and eject the disposable clipper blade
- Easy to clean and disinfect

CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.

DEPARTMENT

Wards/theatre complex



Universal blade



Designed for most hair and cutting needs.

Neuro blade



Designed for coarse or dense hair in areas such as scalp or chest.

Battery charger



Easily mounted on the wall.

MICARE Nitrile Examination Gloves

CODE: NITR-GLOVE-S/M/L/XL

There is an undeniable difference when you put on MICARE Nitrile Examination Gloves. The outstanding comfort and protection gives you the confidence you need to perform your daily tasks.



QUANTITY
200 pcs/pack



Aesthetically Presentable
Unlike other gloves, our unique glove packing method eliminates glove creases – giving you an ironed finish.

Finger Tip Textured
All our nitrile gloves are finger textured for that added handling precision.

Proven Protection
Gloves are clinically tested to handle against a wide range of chemicals which also comply with food safety standards.

Hygienic
Our fully automated glove production process eliminates any possible human contamination from production down to packing.

Comfortable & Stretchy
The ability of our gloves to not only expand, but to also retract according to your hand movements rewards you with a free range motion.


Extended Version
Longer cuff available for that extra protection against splashes over the wrists and arms.

Beaded Cuff
Increased strength & ease in donning.




How Goodpac™ technology dispenses one glove at a time.

With an alternating layering method, each top layer glove cuff will always completely cover the fingers of its immediate bottom layer glove. Hence, exposing only the top most layer for the wearer to pull. Therefore preventing accidental pulling of the bottom layered gloves. This results in zero wastage. Moreover, a neatly packed glove box will prevent creased gloves therefore speeding up the process of dispensing and donning.



When packing, no overlap occurs as gloves are neatly stacked to fit the box



■ 1st layer ■ 2nd layer ■ 3rd layer ■ 4th layer

MICARE Glove Technical Specification

Standards of Conformance

Product conforms to the following test standards:

ASTM D6319	Standard Specification for Nitrile Examination Gloves
ASTM D412	Tensile Properties of Vulcanised Rubbers
ASTM D573	Accelerated Aging for Latex (at 70°C for 166+-2hr)
ASTM D5151	Standard Test for Detection of Holes in Medical Gloves
ASTM D6124	Standard Test Method for Residual Powder on Medical Gloves
ISO 2859-1	Sampling Procedures for Inspection by Attributes
	Part 1:
	Sampling Schemes Indexed by Acceptance Quality Limit (AQL) for Lot-by-lot Inspection

PERFORMANCE REQUIREMENTS

CHARACTERISTIC	SAMPLING PLAN	INSPECTION LEVEL	AQL
Watertight Test	Multiple Normal	G1	1.5
Visual Test	Single Normal	S4	Major 2.5 Minor 4.0
Physical Dimension	Single Normal	S2	AQL 4.0
Physical Properties	Single Reduced	S2	AQL 4.0
Particulate Residue	N/A	N=5	<2mg glove

PHYSICAL DIMENSIONS

SIZE	EX SMALL	SMALL	MEDIUM	LARGE	EX LARGE
Palm width \pm 4 mm	76	86	98	107	115
Glove length (mm)	Min 230				
Thickness \pm 0.02 mm Palm (center of palm) Finger (13 \pm 3 mm from tip)	0.07 0.09				

PHYSICAL PROPERTIES

TEST	BEFORE AGING	AFTER AGING
Tensile Strength (MPa)	Min 18.0	Min 16.0
Elongation (%)	Min 500	Min 400

Nitrile Examination Gloves

CODE: PFSN-FTP-S/M/L/XL

These gloves are suitable for all medical and non-medical applications. They are ambidextrous, i.e. they fit on either left or right hand and are made of nitrile rubber (NBR). Usually 240 mm long, their individual characteristics are:

- Purple blue colour
- Thickness: 3 g
- Powder-free
- Latex-free
- Beaded cuff

These specifications differentiate them in performance, functionality and cost in meeting individual needs and requirements.



QUANTITY
100 pcs/pack



Finger tip textured
for added handling
precision.



Comfortable and stretchy:
ability of the gloves to not
only expand, but to also
retract according to your
hand movement.

Extended version:
longer cuff available
for extra protection.



Protective Coverall

CODE: COVERALL-HJ

The disposable protective coverall is a critical item to provide a barrier and protective work of the medical personnel. It's suitable for operating rooms, medical clinics, hospital wards, inspection rooms, laboratories, ICU & CDC sites for important isolation of virus damage.



QUANTITY

40pcs/pack



Isolation Gown

CODE: GOWN-ISO1000-Y

Front chest, sleeves and back panel are manufactured from SPP (spun bonded polypropylene) material. Elastic wrist cuffs manufactured from 100% polyester. Tie closures at the neck and waist.

ADVANTAGES

- Dust-proof
- Isolate virus
- Disposable
- Environmentally-friendly and hygienic
- Reduce risk by providing level 2 barrier protection during routine care and procedures.
- Can be easily removed and disposed of to lower risk of cross-contamination.



QUANTITY

10pcs/pack



Linen Savers

510 mm x 650 mm x 5-PLY

CODE: LINSAB-65X51X5

Disposable linen saver designed to protect linen from high levels of fluids, providing a waterproof barrier to protect your linen and bedding from contamination. Consists of highly absorbent tissue, sandwiched between a waterproof PE backing sheet and a soft, non-woven fabric.

FEATURES:

- Used for bed protection as well as keeping patients dry, thus avoiding skin irritation
- The highly absorbent pad is fast acting and retains liquid more effectively
- Disposable, single use only

SPECIFICATIONS

- 5-ply
- Waterproof
- Ideal for home or hospital-based care
- Economical



DIMENSIONS

Available in one size only: 510 mm x 650 mm



QUANTITY

10pcs/pack





Surgical Scrub Brush with Nail Cleaner Sterile

CODE: SCRUBB-WNC-MI

Single use for scrubbing of hands and nails in hospital theatre. Soft bristles and sponge with nail cleaner included. Consists of latex-free, dry formulation and easy to open packaging.

ADVANTAGES

- Soft bristles and sponge
- Nail cleaner
- Latex-free
- Dry formulation
- Easy to open packaging



FUNCTION

Single use for scrubbing of hands and nails in hospital theatre.



QUANTITY

300 pcs/pack



Manual Alcohol/Soap Dispenser including Drip Tray & Pump

CODE: DISP-MAN-00

This quality designed ABS plastic dispenser can accommodate a 800ml liquid sachet or 800ml refillable plastic container and uses an easy push mechanism for the dispensing of liquid sanitiser or soap. This item comes equipped with a drip tray and pump.



FUNCTION

Manual dispensing of one 800 ml liquid bag containing alcohol or soap solution.



DIMENSIONS

H 375 mm x W 125 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Ideally for the home or office.
Can be used at entrances or any high traffic area.

- Contents of fluid solution plastic bag should not exceed 800 ml (range 600 ml - 800 ml)
- Fluid solution plastic bag should be made from durable material, to withstand handling as well as the effects of the fluid content
- Fluid solution plastic bag outlet should be positioned at the bottom edge of the bag
- Valve inlet inner diameter of fluid solution plastic bag: 6.8 mm



Outer dimensions of the fluid solution plastic bag should be no larger than 196 mm x 246 mm



Also available: Refillable plastic container



Drip Tray



Pump



Manual Alcohol/Soap Dispenser Only

CODE: DISP-MAN-01

This quality designed ABS plastic dispenser can accommodate a 800ml liquid sachet or 800ml refillable plastic container and uses an easy push mechanism for the dispensing of liquid sanitiser or soap. This item comes equipped with a pump.

- Contents of fluid solution plastic bag should not exceed 800 ml (range 600 ml - 800 ml)
- Fluid solution plastic bag should be made from durable material, to withstand handling as well as the effects of the fluid content
- Fluid solution plastic bag outlet should be positioned at the bottom edge of the bag
- Valve inlet inner diameter of fluid solution plastic bag: 6.8 mm



FUNCTION

Manual dispensing of one 800 ml liquid bag containing alcohol or soap solution.



DIMENSIONS

H 375 mm x W 125 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Ideally for the home or office. Can be used at entrances or any high traffic area.



Outer dimensions of the fluid solution plastic bag should be no larger than 196 mm x 246 mm



Also available: Refillable plastic container

Manual Alcohol/Soap Dispenser with Refill Container

CODE: DISP-MAN-05

This quality designed ABS plastic dispenser can accommodate a 800ml liquid sachet or 800ml refillable plastic container and uses an easy push mechanism for the dispensing of liquid sanitiser or soap. This item comes equipped with a refillable container and pump.



DIMENSIONS

H 375 mm x W 125 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Ideally for the home or office. Can be used at entrances or any high traffic area.



Refill Container

Pump included



Alcohol/Soap Touchless Dispenser

CODE: DISPENSER-03

The dispenser consists of an ABS plastic dispenser case and control box. This touchless automatic dispenser accommodates one 800 ml liquid bag containing alcohol or soap solution. Available liquid solutions: Anti-bacterial hand soap, hand sanitiser. Built-in electronic infrared (IR) sensor activates the dispensing cycle. The dispenser is powered by 6 x AA batteries (not included in price). The battery lifespan is approximately 10 000 cycles. A toggle switch can be adjusted to dispense either one (0.4 ml), two or three drops of liquid.

ADVANTAGES

- Convenient, reliable hands-free operation.
- Reduces operating costs due to fixed liquid dispensing.
- Time saver.
- Eliminates risk of germ transmission between users.
- Easy liquid bag replacement.
- Lockable for added tamper-resistance.



FUNCTION

Touchless automatic dispenser accommodates one 800 ml liquid bag containing alcohol or soap solution. The dispenser can be wall mounted or mounted to a utility pole.



DIMENSIONS

H 260 mm x W 125 mm x D 135 mm



CLEANING

To be cleaned with a disinfectant agent, NOT a chlorine-based chemical.



DEPARTMENT

Suitable for use in homes, kitchens, restaurants, hotels, offices, beauty salons, schools and public areas.

Liquid Hand Sanitiser 800 ml Sachet

CODE: PRIMESTPAR

Dispensing sachet for hands-free or manual sanitiser, containing 80% ethanol-based hand sanitising liquid.



Dispenser Pump

CODE: DISP-MAN-03

Accessory to hands-free and manual dispenser



Dispenser Refill Container

CODE: DISP-MAN-01

Accessory to manual dispenser



Dispenser Drip Tray

CODE: DISP-MAN-02

Accessory to hands-free and manual dispenser

T: +27 (0)21 851 8484

E: sales@medicalinnovations.co.za

W: medicalinnovations.co.za

7 Cyclonite St, The Interchange
Somerset West, 7130
Western Cape, South Africa

The information and dates supplied in this catalogue, are based on the information believed to be reliable, however are not binding and inaccuracies may occur. Medical Innovations, reserves the right to change specification and design at any time without notice.

© 2017 by South Africa. All rights reserved under the Copyright Act of 1987 with the amendments. Reproduction of any of the content, full or partial, is allowed with prior written approval of South Africa. No reasonable requests will be refused provided the condition as laid out are met.